

TSIII SOI Fixture 06.2016



- Submerged type implant with an internal hex 11° tapered connection structure
- Super hydrophilic surface with superior blood wettability, coated with K material
- Super hydrophilic surface inducing fast blood clot formation
- Tapered body design for excellent initial stability
- Excellent initial stability in soft bone with smaller threads in the upper section
- Superior self-threading effect with corkscrew thread
- Excellent initial stability needed for immediate loading even in soft bone
- Recommended placement torque : $\leq 40\text{Ncm}$
- ※ Fixtures with a diameter of 4.5mm or greater are recommended for the posterior region.

Narrow


- Used in tight spaces (narrow ridge)
- Easy angle compensation in anterior region
- Compatible with existing mini upper parts (not compatible with cover screw, mount, or lab analog)

NoMount Fixture order code

: fixture product code (ex : TS3S4010A)

Pre-Mounted Fixture (fixture + mount + cover screw) order code


: B + fixture product code (ex : B**TS3S4010A**)




D \varnothing 4.5
Hex 2.5



D \varnothing 5.0
Hex 2.5




TS3S4507A TS3S4508A TS3S4510A TS3S4511A TS3S4513A



TS3S5007A TS3S5008A TS3S5010A TS3S5011A TS3S5013A


01.2019

D \varnothing 3.0
Hex 2.1
Narrow




TS3M3008A TS3M3010A TS3M3011A TS3M3013A

D \varnothing 3.5
Hex 2.1



TS3M3508A TS3M3510A TS3M3511A TS3M3513A

D \varnothing 4.0
Hex 2.5



TS3S4007A TS3S4008A TS3S4010A TS3S4011A TS3S4013A

Nominal diameter may differ from actual diameter.