


Transfer Abutment 01.2013

- Abutment for producing cement-retained/combination prosthesis
- Fixture level impression
- Abutment level impression possible by rigid impression coping (except Ø4.0)
- Tightened with 1.2 hex driver
- Recommended tightening torque : 20Ncm (mini), 30Ncm (regular)
- Packing unit : Abutment + EbonyGold screw



Abutment + EbonyGold screw order code
: Product code + **WH** (ex : GSTA4621**WH**)

D Ø4.0



EbonyGold screw
: GSABSM

Abutment level impression not available

		H \ G/H	1.0	2.0	3.0	4.0
Hex	5.5					
	7.0		GSTA4612	GSTA4622	GSTA4632	GSTA4642
Non-Hex	5.5					
	7.0		GSTA4612N	GSTA4622N	GSTA4632N	GSTA4642N
		H \ G/H	5.0	6.0	7.0	
Hex	5.5					
	7.0		GSTA4652	GSTA4662	GSTA4672	
Non-Hex	5.5					
	7.0		GSTA4652N	GSTA4662N	GSTA4672N	

D Ø4.5



EbonyGold screw
: GSABSM

		H \ G/H	1.0	2.0	3.0	4.0
Hex	5.5					
	7.0		GSTA4611	GSTA4621	GSTA4631	GSTA4641
Non-Hex	5.5					
	7.0		GSTA4611N	GSTA4621N	GSTA4631N	GSTA4641N
		H \ G/H	5.0	6.0	7.0	
Hex	5.5					
	7.0		GSTA4651	GSTA4661	GSTA4671	
Non-Hex	5.5					
	7.0		GSTA4651N	GSTA4661N	GSTA4671N	

D Ø4.5



EbonyGold screw
: GSABSS

		H \ G/H	1.0	2.0	3.0	4.0
Hex	5.5					
	7.0		GSTAS4611	GSTAS4621	GSTAS4631	GSTAS4641
Non-Hex	5.5					
	7.0		GSTAS4611N	GSTAS4621N	GSTAS4631N	GSTAS4641N
		H \ G/H	5.0	6.0	7.0	
Hex	5.5					
	7.0		GSTAS4651	GSTAS4661	GSTAS4671	
Non-Hex	5.5					
	7.0		GSTAS4651N	GSTAS4661N	GSTAS4671N	

Transfer Abutment 01.2013

D Ø5.0



EbonyGold screw
: GSABSS

		H \ G/H	1.0	2.0	3.0	4.0
Hex	4.0					
	5.5					
	7.0					
Non-Hex	4.0					
	5.5					
	7.0					

D Ø6.0



EbonyGold screw
: GSABSS

		H \ G/H	5.0	6.0	7.0
Hex	4.0				
	5.5				
	7.0				
Non-Hex	4.0				
	5.5				
	7.0				

		H \ G/H	5.0	6.0	7.0
Hex	4.0				
	5.5				
	7.0				
Non-Hex	4.0				
	5.5				
	7.0				

D Ø7.0



EbonyGold screw
: GSABSS

		H \ G/H	1.0	2.0	3.0	4.0
Hex	4.0					
	5.5					
	7.0					
Non-Hex	4.0					
	5.5					

		H \ G/H	5.0	6.0	7.0
Hex	4.0				
	5.5				
	7.0				
Non-Hex	4.0				
	5.5				

D Ø6.0



EbonyGold screw
: GSABSS

		H \ G/H	1.0	2.0	3.0	4.0
Hex	4.0					
	5.5					
	7.0					
Non-Hex	4.0					
	5.5					
	7.0					