



# Digital Prosthetics Scan Healing Abutment

## TS Scan Healing Abutment NEW 2021

- Healing Abutment with scan body function
- Different types are differentiated by the marker shape at the top (See the type classification table)
- Abutment level impression
- Driver-integrated carrier is used for transfer to the mouth and ease of connection
- Different screws are used for different lengths (mixed use not possible)
- 1.2 hex driver / driver-integrated carrier
- Packing unit : scan healing abutment body + Ti screw

Scan healing abutment body + Ti screw order code  
: product code + TH (ex : TSSHA505RTH)

- M** Mini
- R** Regular




D \ H	4.0	5.0	7.0	9.0
<b>Ø4.0</b>	 TSSHA404M	 TSSHA405M	 TSSHA407M	 TSSHA409M

D \ H	4.0	5.0	7.0	9.0
<b>Ø4.5</b>	 TSSHA454R	 TSSHA455R	 TSSHA457R	 TSSHA459R
<b>Ø5.0</b>	 TSSHA504R	 TSSHA505R	 TSSHA507R	 TSSHA509R
<b>Ø6.0</b>	 TSSHA604R	 TSSHA605R	 TSSHA607R	 TSSHA609R

## TS Scan Healing Abutment Carrier NEW 2021

- Used for transferring the Scan Healing Abutment into the mouth and for connection
- Selected according to the body diameter for use
- Material: PEEK + TrimRite

- M** Mini
- R** Regular



Type \ Scan Healing(D)	Ø4.0	Ø4.5	Ø5.0	Ø6.0
<b>Short</b>	 TSSHAC400	 TSSHAC450	 TSSHAC500	 TSSHAC600
<b>Long</b>	 TSSHAC400L	 TSSHAC450L	 TSSHAC500L	 TSSHAC600L